herpotherm®

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herpotherm®

simple | quick | well-tolerated

A medical device for the symptomatic treatment of cold sores

Please read the entire instructions for use carefully before the first application and follow the instructions and warnings.

Intended use

herpotherm® is a medical device for the local treatment of symptoms that accompany cold sores (herpes simplex labialis) such as itching, burning, tingling or tightness.

herpotherm® can normally prevent the development of cold sore blisters if applied early enough. If applied later, the severity of symptoms may be lessened, or they may fade sooner. The first signs of an outbreak (so-called prodromal symptoms) are generally itching, burning, tingling or tightness.

How it works

herpotherm® works solely by concentrated heat and thus completely free from chemicals.

Hence, herpotherm® is also suitable for pregnant women, allergy sufferers and children.


The application of concentrated heat (local thermotherapy) is a physical mode of action based on a brief, concentrated application of heat to a small, limited area of skin. Even a localised pulse of heat may be sufficient to relieve symptoms such as itching, burning, tingling or tightness and, if applied early enough to normally prevent the formation of cold sore blisters.

Application

herpotherm® is very simple to use:

- Remove the protective cap from the device.
- Place the ceramic contact plate of the device onto the affected skin area and gently press it against the skin.
- Press the button and the ceramic contact plate will heat up to a temperature of around 51°C. Once this temperature is reached, a LED will light up. The LED will go off once the device automatically shuts down after 3 seconds of use. At the start of the heating-up phase and at the end of treatment, you will hear an additional acoustic signal.
- Place the protective cap on the device.

For persistent or recurrent symptoms, multiple treatment of the same skin site is possible. Please leave a treatment interval of at least 2 minutes between each application. Do not exceed a maximum number of 5 applications per hour on the same skin area.

Directions for use

herpotherm® is indicated for self-administration by adults and children aged 12 years and older.

In children under 12 years of age or in people with significant impairment of hearing, vision and/or movement, the device must only be used with the assistance of an adult.

Do not use the device in children under 3 years of age. Due to the cognitive development status respectively the development of skin tissue the treatment with the device is recommended for children from the age of 7.

To prevent virus transmission and for reasons of hygiene, herpotherm® should only ever be used by one person.


Cold sores are a viral disease caused by the herpes simplex virus. It is transmitted by droplet infection, by direct contact or as a smear infection.

There is no permanent cure for infection with the herpes simplex virus.

Cleaning information:

The ceramic contact plate of the device can be easily cleaned with a standard alcohol disinfectant wipe; all other external surfaces, for example, with a clean cloth.

Battery replacement:

To change the battery, the protective cap and battery cap must be removed and refitted after inserting the battery (see illustration). When attaching the battery cap, make sure that the „Herpotherm®“ logo and the flattened edge of the battery cover have the same orientation so that the housing can be closed properly.

When the battery needs changing, the LED will flash 3 times and the acoustic signal will sound 3 times.

Do not use rechargeable batteries.

Warnings

- herpotherm® must not be used in or on
 - eyes and eyelids
 - mucous membranes
 - open wounds, damaged skin (e.g. open blisters)
 - scalded or burnt skin
 - frostbite
 - known increased allgesia in small fibre neuropathy, allodynia, hyperalgesia
 - tumours and malignant changes of the skin
 - people with known limited sensitivity to pain and temperature.

Medical advice should be sought in the following cases:

- persistent symptoms while using herpotherm®.
- worsening symptoms while using herpotherm®. Until these have been clarified, you should stop using the device.

Furthermore, the device must not be used:

- immediately after application of ointments, creams, lotions etc. on the affected site
- on moist skin.

In individual cases, the treatment temperature may cause some discomfort. On sensitive skin, a temporary reddening of the skin may occur after use. In very rare cases, skin irritation may occur in users with a sensitive skin type.

The device indicates failure by flashing the LED 5 times and sounding the acoustic signal 5 times. This may occur if usage is incorrect, e.g. if the user triggers and operates the device in the air instead of on the skin. If this error message occurs repeatedly despite proper use, please contact the manufacturer or distributor.

Do not use the device in case of damage, obvious malfunction or with open housing. Do not use batteries that are damaged or moist. Protect the device from moisture.

Permissible environmental conditions

Store the device protected from direct sunlight, in a dry and clean place.

Please observe the following permissible environmental conditions.

For operation	- Environmental temperature 5 °C to 40 °C - Relative humidity 15 % to 90 %, non-condensing - Air pressure of 700 hPa to 1060 hPa
For transport and storage without packaging*	-25 °C to 5 °C -5 °C to 35 °C for a relative humidity of up to 90 %, non-condensing -35 °C to 70 °C for a water vapour pressure of up to 50 hPa - Air pressure of 700 hPa to 1060 hPa
For transport and storage with packaging*	- Environmental temperature 5 °C to 40 °C - Relative humidity 5 % to 85 % - Air pressure of 700 hPa to 1060 hPa
* Please note that the manufacturer of the batteries may have stated different storage conditions.	

Do not use the battery beyond the expiry date which is stated on the battery. In order to avoid damage to herpotherm® due to leaking battery acid, the battery should be taken out during prolonged periods of non-use (inactivity). Please observe the handling instructions on the batteries you use. Dispose of used batteries according to local regulations for disposal of batteries.

Please do not dispose of the device in normal household waste. Dispose of the device at the end of its expected service life of 8 years in accordance with local regulations for the disposal of waste electrical equipment.

Manufacturer:
© mibetec GmbH
Münchener Straße 15
06796 Brehna
Germany

www.herpotherm.com

Distributor:

Mibe Pharma UK Ltd
4 Coleman Street
6th floor
London
United Kingdom, EC2R 5AR
Tel.: +44(0)2038707504
e-mail: customerservice.mibeUK@dermapharm.com
www.myherpotherm.co.uk

Please contact the manufacturer or distributor if you have any questions about the device or how it works. Please report any serious incident that has occurred in relation to herpotherm® to the manufacturer and the competent authority in your country.

Date of revision of the text: 2022-06

CE 0123

	Please read the instructions for use	REF	Catalogue number
MD	Medical Device	LOT	See device for batch number
	Protect the device from direct sunlight and store in a dry place/protect from moisture.		Do not dispose of with household waste/dispose according to local regulations for disposal of waste electrical equipment
 40 °C 5 °C	Permissible environmental temperature for transport and storage with packaging	 1060hPa 700hPa	Technical data: Battery type/-voltage CR2 3 V
 85 % 5 %	Permissible relative humidity for transport and storage with packaging		The ceramic contact surface is a type BF applied part in accordance with DIN EN 60601-1

This medical device complies with the EMC standard IEC 60601-1-2:2014.

Nevertheless, special precautions should be taken:

herpotherm® is intended for use in home healthcare environment. Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation.

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm to any part of herpotherm®, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

According to the specifications of IEC 60601-1-2:2014, the technical description can be requested in digital form and English or German language from mibetec GmbH (address can be found in the instructions for use).

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